MA2701

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Schaukel-Massagesessel

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Instruction Manual

Rocking Massage Chair

Please read this manual carefully before use
and keep it safe for future reference.

Common problem handling methods

If you encounter any problems while using this product, please refer to the following tips and possible solutions.

If the product remains malfunctioning or unusable, please turn off the power and unplug it, and contact our service center.

Situation	Possible causes/solutions
Not working	● Please ensure the switch is turned on. ● Please check that the power plug is properly inserted.
The massage suddenly stopped during the process.	● The button on the manual controller was accidentally pressed. ● The 15-minute timer set by the automatic timer has elapsed. ● If the overheat protection device is activated, please turn off the main power and wait for the massager to cool down. ● If the problem persists, please do not continue to use it and contact the company's service center.

Technical Parameters

Product Name: Rocking Massage Chair

Product Model: MA2701

Rated Power: 36W

Rated Voltage: 12V

The company reserves the right to improve the product and instruction manual. Changes will be made without prior notice. The actual product appearance shall prevail.

This booklet has been carefully checked, but there may be printing errors or misunderstandings in the content. We appreciate your understanding and reserve the right of interpretation.

♻ Improper disposal of waste can cause environmental pollution.

Certificate of Conformity

Product Name: Rocking Massage Chair

Model: MA2701

PASS

QC: Inspection 08

PASS

Production date:

Precautions

When using this massage product, especially in the presence of children, the following safety precautions must be observed. To avoid the risk of electric shock, please read all instructions carefully before use.


- Do not use in humid environments. Do not place this product in water or other liquids.
- Do not use needles or other metal objects to secure this product. Do not use unrecommended accessories or attachments. Do not use this device if the casing or power cord is damaged.
- If the product gets wet, do not use it yourself. Contact the manufacturer immediately for after-sales service.
- Do not use this device if you are intoxicated or feeling unwell. Do not use this product within one hour after meals. If the massage area is injured or you have recently undergone surgery, please use it under the guidance of a doctor.
- If the machine malfunctions or does not work properly, stop using it immediately and take it to the store where it was purchased or an authorized service center for repair. Do not attempt to repair it yourself.
- Do not use for extended periods; do not exceed 30 minutes.
- This product is not recommended for use by persons with physical disabilities, sensory or mental impairments, or lack of relevant experience or knowledge (including children) unless supervised or instructed by a person responsible for their safety.
- Children should be supervised to ensure they do not play with this product.
- This product has a heating surface; caution must be exercised when using it if you are not sensitive to heat.
- This product is powered by the provided power adapter. If damaged, a replacement power adapter of the appropriate specification must be purchased from the manufacturer, its service department, or a similar institution.

- Do not pull the connecting cord directly when plugging or unplugging the power cord.


WARNING:
Consult a doctor before use if you are undergoing treatment or have any health conditions, or if you have any of the following conditions: heart disease patients, malignant tumor patients, pregnant women, women during menstruation, people with spinal deformities or abnormalities, people with a history of treatment using the product on affected areas, people with osteoporosis, or patients with implanted pacemakers.

Connection diagram

Connect the DC plug of the adapter to the connector on the right side of the massage rocking chair, as shown in the picture, and connect the other end to the power socket.

Connection diagram

⏻ : One-button power-on, defaulting to full-range back and waist massage with hip heating.

Full Back: Turns the full-range back and waist massage on or off. When the massager is on the upper or lower back, press this button to switch to full-range back and waist massage.

Upper Back: Press once to massage the upper back area, press again to turn off the massage, and repeat.

Lower Back: Press once to massage the lower back (waist) area, press again to turn off the massage, and repeat.

Mode: Neck massage mechanism forward/reverse adjustment. Default is forward rotation upon startup. Press once to reverse rotation, press again to turn off the neck massage, and repeat.

Heating: Default is hip heating massage upon startup. Press once to turn off heating massage, and repeat.

Intensity: Massage mechanism speed adjustment. Default is medium speed upon startup. Press once for low speed, press again for high speed, and repeat.

Vibration: Press once to activate hip vibration massage. Default is low vibration speed. Press again for medium vibration speed, press again for high vibration speed, press again to turn off vibration massage, and repeat.